Warum TESTA

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Why TESTA

TESTA Analytical Solutions is fully focused on Polymer Characterization and Particle Sizing Technologies, offering the most advanced and refined equipment for the most challenging research, development and QC Tasks. TESTA Analytical Solution's experts are known for their ability to help users solve the most intricate problems and applications, worldwide. Our know-how is the fundament of our success.

RELATED PRODUCTS

AYE dn/dc Refractometer

AYE Differential RI & Viscometer Combo

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
AYE Differential Refractive Index Detector

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Why a DRI

A Differential Refractive Index Detector (DRI) is the most flexible and widespread detector in GPC/SEC Chromatography. It provides high sensitivity, a broad linear range and is capable to detect (almost) any dissolved compound. Due to this, it has been the detector of choice for many years. It is utilized as the sole detector in a simple GPC/SEC system or as the source of concentration data in a more sophisticated setup including a Viscosimeter and/or a MALS Detector. Its applicability is nearly infinite, as such the DRI is undoubtedly the most useful detector for most polymer characterization applications.

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Why the AYE DRI

GPC/SEC Chromatography with MALS Detector is a very powerful tool capable of delivering absolute values of the compound under investigation. To achieve this, it is absolutely necessary to obtain concentration data along with the molar mass. In this scenario, a common DRI will lack the capability to deliver data measured at the same wavelength as the laser used in the MALS Detector. This fact will introduce a non-neglectable error in all calculations performed, thus questioning the obtained results. The AYE DRI operates at the very same wavelength your MALS Detector does, thus improving reliability of your results.