


Versatile UV Spectrophotometer for batch and flow applications

Uniqsis reports that demand for its **Flow-UV inline UV-Visible spectrophotometer** for **biological, biotech, drug discovery and industrial applications** grew to record levels during 2022.



Compact in design, the Flow-UV from Uniqsis sets a benchmark for real time monitoring of continuous flow and batch applications. Unlike many UV-visible spectrophotometers the high-resolution Flow-UV CCD array detector does not require calibration or routine servicing. The Xenon flash lamp source used in the Flow-UV has a lifetime of up to 10 years.


Steve Evans,

Commercial Director for Uniqsis commented “During the last 12 months, one biotech customer ordered 6 Flow-UV units for nucleic acid purity assessment. In a rapidly growing number of pharma labs the unit is being used for high throughput screening of UV-active drug substances and real time determination of the steady state concentration of products being produced in drug discovery reaction optimisation studies. Researchers are also using the Flow-UV for monitoring algal growth, looking at variations of different metal ion complex solutions in rock cores and optimising microporous membrane production to name but a few of the wide array of applications this versatile spectrophotometer was selected for”.



He added

“The Flow-UV uses fibre optic waveguides that connect the flow cell to the detector, allowing the flow cell to be positioned virtually anywhere in the flow path. To assure linearity of response, up to 5 wavelengths can be selected to monitor the progress of a reaction or monitor concentration changes in a wide variety of media. The Flow UV is supplied with 2 software applications for both flow (dynamic) and batch processes. Customers really appreciate the intuitive method set up and simple data acquisition interface.”



For further information

about the Flow-UV please visit www.uniqsis.com/paProductsDetail.aspx?ID=Flow%20UV%20ACC or contact Uniqsis on +44-845-864-7747/ info@uniqsis.com.

Since 2007, Uniqsis has specialised in the design and supply of mesoscale continuous flow chemistry systems for a wide range of applications in chemical and pharmaceutical research. The company's aim is to make flow chemistry easily accessible to both novices and experienced users.

Further Information:

Uniqsis Ltd

Foxes Bridge Farm
Unit 1, Lime Tree Barn
Foxes Bridge Farm
Royston Lane
Comberton
Cambridgeshire CB23 7EE
UK

tel: +44-1223 942004

email: info@uniqsis.com

web: www.uniqsis.com